



# IRiS

## Colour Selector

### Instruction Manual



#### **WARNING:**

Choking Hazard – Small Parts  
Not suitable for children under 3  
years.



3+ Years

**experia**  
innovation in every sense

**Congratulations** on the purchase of your IRiS Colour Selector!

The IRiS Colour Selector has been designed specifically for people with special needs and should, if used in accordance with these guidelines, offer years of switching learning and fun!

The IRiS Colour Selector offers endless opportunities to develop numerous skills including colour recognition, cause and effect, hand to eye co-ordination, gross and fine motor skills and many more.

The IRiS Colour Selector is an IRiS Talker and can control any number of IRiS Listeners.

Put simply, it is operated by touching any of the switches.

**Please ensure you read this Instruction Manual in full before operating the unit.**

### Features

- IRiS Wirefree Compatible
- 8 Brightly Coloured Switches
- 3.5mm Switch Input Sockets
- Capacitive Technology for Touch Sensitive Control
- 4 x AA Batteries
- Easy to use
- 2 Year Manufacturers Guarantee (Collect and Return)

### Getting Started

1. Simply touch any switch on the IRiS Colour Selector – the unit will control the products it is paired to.
2. The unit features a power saving mode to prolong battery life. Touch any of the following colours initially to wake the product: Red, yellow, Green, Blue.
3. To pair the IRiS Colour Selector to a product please see the 'Pairing to IRiS Listeners' section in this manual.

### Installing Batteries

1. Remove the battery cradle from the unit by lifting each edge as shown on the compartment.
2. With the compartment removed, insert batteries as shown within the compartment.
3. Push the battery compartment into place ensuring your hands **DO NOT** touch the front of the unit.
4. Wait 10 seconds and the unit will be ready to use.

**Please note:** this unit uses capacitive technology to make the surface on the front of the unit touch sensitive. When the unit is initially powered it calibrates itself ready for use. If you have a hand or finger on the front of the unit on start up, it unit will **NOT** calibrate correctly. If this does happen, please remove battery compartment, wait 10 seconds then replace batteries ensuring nothing is touching the front of the unit.

### Using your own switches

Connect your own switches by plugging them into the sockets at the rear of the unit. The sockets are located in the same order as the colours on the top of the unit. 8 x 3.5mm jack adaptors are integral enabling you to use your own switch.

### Pairing to IRiS Listeners

This procedure is split into 2 sections, one for those who do not have an IRiS Master Pairer and for those who do.

#### *Without a IRiS Master Pairer*

Every IRiS Listener has a button on it called a pairing button. To get the IRiS Colour Selector to work with a product, first decide which product you want to control then find the pairing button on that product and press it. You then need to press any switch on IRiS Colour Selector within 5 seconds. The IRiS Listener will now change colour every time you press a switch if correctly paired. The unit will now work with the IRiS Colour Selector.

#### *With a IRiS Master Pairer*

To get the IRiS Colour Selector to work with a product, first decide which unit you want to control and find its pairing button on the IRiS Master Pairer and press it. You then need to press any switch on IRiS Colour Selector within 5 seconds. The IRiS Listener will now change colour every time you press a switch if correctly paired. The unit will now work with the IRiS Colour Selector.

### Cleaning and Maintenance

The unit should be cleaned with a damp soft cloth. A non abrasive cleaning agent may be used for stubborn stains. Ensure no liquid enters the unit.

The unit has been designed for use by people with special needs, however this does NOT mean it is unbreakable. You must ensure the unit is not dropped, thrown, stood upon, immersed in water, etc.

### Troubleshooting

Q. The unit doesn't work.

A2. Replace the batteries

A3. If the unit does not operate correctly after completing the above, please contact us.

Q. The unit works sporadically.

A1. Replace batteries making sure you do NOT touch the front part of the unit.

A2. If the unit does not operate correctly after completing the above, please contact us.

Q. The unit won't control any product.

A1. When a switch is pressed on the unit does a blue LED illuminate? If not,

A2. Replace the batteries

A3. Is the unit paired to an IRiS Listener? If so,

A4. Try pairing the unit to another IRiS Listener, does this work? If not, please contact us.

### Contact Details

For USA please email [support@experia-usa.com](mailto:support@experia-usa.com).

For UK and rest of the world email [support@experia-innovations.co.uk](mailto:support@experia-innovations.co.uk)

Telephone: 0044 (0) 845 644 0977.

Please note: email is our support teams preferred communication and this takes precedence over telephone calls.

### QUESTIONS & IDEAS?

If you have an idea on how to improve this product or want to see your product idea realised, please email: [ideas@experia-innovations.co.uk](mailto:ideas@experia-innovations.co.uk).

All ideas are presented to our OT, technical and finance department and if all procedures are passed we will contact you to discuss your idea further.

We can't promise to develop every product but if it's a great idea, we'll try!

### **Manufactured proudly in the UK by:**

Experia

Acorn Phase 3, High Street, Grimethorpe, Barnsley, S72 7BD, South Yorkshire

